




## 谷胱甘肽 S-转移酶

### PRINCIPLE:



Abbreviations:

G-SH = Glutathione, Reduced Form

CDNB = 1-Chloro-2,4-Dinitrobenzene

G-SDNB Conjugate = Glutathione-2,4-Dinitrobenzene

### CONDITIONS:

T = 25°C, pH = 6.5, A340nm, Light path = 1 cm

### METHOD:

Continuous Spectrophotometric Rate Determination

### REAGENTS:

A. 100 mM Potassium Phosphate Buffer with 1.0 mM Ethylenediaminetetraacetic Acid, pH 6.5 at 25°C (Prepare 100 ml in deionized water using Potassium Phosphate, Monobasic, Anhydrous, and Ethylenediaminetetraacetic Acid, Tetrasodium Salt: Hydrate, Adjust to pH 6.5 at 25°C with 1M KOH)

B. 75 mM Glutathione, Reduced Solution (G-SH)  
(Prepare 10 ml in Reagent A using Glutathione, Free Acid, Reduced Form, PREPARE FRESH.)

C. 30 mM 1-Chloro-2,4-Dinitrobenzene Solution (CDNB)  
(Prepare 1 ml in Reagent D using 1-Chloro-2,4-Dinitrobenzene, PREPARE FRESH.)

D. 95% Ethanol (Nondenatured)

E. Glutathione-S-Transferase Enzyme Solution

(Immediately before use, prepare a solution containing 0.075 - 0.15 unit/ml of Glutathione S-Transferase in cold Reagent A.)



## PROCEDURE:

Pipette (in milliliters) the following reagents into suitable cuvettes:

|                                                                                                                      | Test | Test |
|----------------------------------------------------------------------------------------------------------------------|------|------|
| Reagent A (Buffer)                                                                                                   | 2.70 | 2.80 |
| Reagent B (G-SH)                                                                                                     | 0.10 | 0.10 |
| Reagent C (CDNB)                                                                                                     | 0.10 | 0.10 |
| Equilibrate to 25° C. Monitor the A340nm until constant, using a suitably thermostatted spectrophotometer. Then add: |      |      |
| Reagent E (Enzyme Solution)                                                                                          | 0.10 | -    |

Immediately mix by inversion and record the increase in A340nm for approximately 5 minutes. Obtain the  $\Delta A_{340\text{nm}}/\text{minute}$  using the maximum linear rate for both the Test and Blank.

## CALCULATIONS:

$$\text{Units/ml enzyme} = \frac{(\Delta A_{340}/\text{min Test} - \Delta A_{340}/\text{min Blank})(3.0)(\text{df})}{(9.6)(0.10)}$$

3.0 = Total volume (in milliliters) of assay

df = Dilution factor

9.6 = Millimolar extinction coefficient of Glutathione-1-Chloro-2,4-Dinitrobenzene conjugate at 340 nm

0.10 = Volume (in milliliter) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

## UNIT DEFINITION:

One unit will conjugate 1.0  $\mu\text{mole}$  of 1-chloro-2,4-dinitrobenzene with reduced glutathione per minute at pH 6.5 at 25°C.

## FINAL ASSAY CONCENTRATION:

In a 3.00 ml reaction mix, the final concentrations are 97 mM potassium phosphate, 0.97 mM ethylenediaminetetraacetic acid, 2.5 mM glutathione, reduced,



上海源叶生物科技有限公司  
Shanghai yuanye Bio-Technology Co., Ltd  
电话: 400-666-5481 传真: 021-55068248  
网址: [www.shyuanye.com](http://www.shyuanye.com)  
邮箱: [3008007435@qq.com](mailto:3008007435@qq.com)

1.0 mM 1-chloro-2,4-dinitrobenzene, 3.2% (v/v) ethanol and 0.0075 - 0.015 unit glutathione-S-transferase.

#### NOTES:

1. This assay is based on the cited references.
2. Where Yuanye Product or Stock numbers are specified, equivalent reagents may be substituted.



源叶