



上海源叶生物科技有限公司
Shanghai yuanye Bio-Technology Co., Ltd
电话: 021-61312973 传真: 021-55068248
网址: www.shyuanye.com
邮箱: shyysw@sina.com

磷酸二酯酶I 检测方法

1 原理

The assay is essentially that of Razell and Khorana (1959) where the reaction velocity is determined by an increase in absorbance at 400 nm resulting from the hydrolysis of p-nitrophenyl thymidine-5'-phosphate. One unit hydrolyzes one micromole of p-nitrophenyl thymidine-5'-phosphate per minute at pH 8.9 and 25°C under the specified conditions.

2 试剂:

- 0.11 M Tris· HCl buffer, pH 8.9, with 0.11 M NaCl and 15 mM MgCl₂ (Tris· Salts buffer)
- 5 mM p-nitrophenyl thymidine-5'-phosphate. Note: The purity of commercial preparations varies somewhat and should be considered in preparing this reagent.

3 酶:

Dissolve at 1 mg/ml in Tris*Salts buffer to obtain a rate of 0.02-0.04 ΔA/minute.

4 操作规程:

Set spectrophotometer at 400 nm and 25°C. Pipette into microcuvettes as



上海源叶生物科技有限公司
Shanghai yuanye Bio-Technology Co., Ltd
电话: 021-61312973 传真: 021-55068248
网址: www.shyuanye.com
邮箱: shyysw@sina.com

follows:

Tris·Salts buffer	0.9 ml
5 mM p-nitrophenyl thymidine-5'-phosphate	0.1 ml

Incubate cuvettes in spectrophotometer for 3-5 minutes to reach temperature equilibrium and establish blank rate, if any. Add 10 microliters of diluted enzyme and record increase in A_{400} for 3-5 minutes. The reaction remains linear until A_{400} reaches about 1.2. Calculate $\Delta A_{400}/\text{minute}$ from initial linear portion of absorbance curve.

5 计算

$$\frac{\text{Units}}{\text{mg}} = \frac{\frac{\Delta A_{400}}{\text{min}}}{16 \times \frac{\text{mg enzyme}}{\text{ml reaction mixture}}}$$

(where 16 is the extinction coefficient of p-nitrophenol under these conditions)