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INULINASE ACTIVITY QUANTIFICATION

Aim:

To check the inulinase activity of the powder or liquid inulinase formulations.

Methods and Materials:

The DNS method was used to determine the enzyme activity in the product, and fructose was used as the substrate.

1. Preparation of the reagents

1.1. Liquid DNS (3,5-Dinitrosalicylic acid) agents preparation:

Solution A-- Dissolved 6.3 g DNS in the water at 45 °C. Add 262 mL 2M NaOH solution, and stir till the solution becoming clear.

Solution B—Dissolved 182 g potassium sodium tartrate in 50 mL water.

Mix Solution A and Solution B together, and then add 5 g crystalline phenol and 5 g NaHSO₃ (no water) till them dissolved totally (yellow color).

The DNS work solution was kept at room temperature for about 7 days to make it stable and then could be used for the following inulinase activity determination. The work solution will be stable for half a year.

1.2. Standard fructose solution (1 mg/mL) preparation.

Dissolved 0.50 g fructose into the distilled water to make the volume to be 500 mL in volumetric flask. The final concentration was 1 mg/mL.

1.3. Preparation 5% inulin solution.

1.3.1. HAc-NaAc (pH4.6, 0.1 M) buffer.

1) 2 M NaAc stock solution: Dissolved 272.16 g NaAc into distilled water and make the final volume was 1000 mL.

2) 2 M HAc stock solution: Dissolved 120.10 g into the distilled water and make the final volume was 1000 mL.

3) Aliquot 24.5 mL 2M NaAc stock solution and 25.5 mL 2M HAc stock solution, and then mix together. Add distilled water to make the final volume as 1000 mL. The 0.1 M HAc-NaAc (pH 4.6) buffer was prepared.

1.3.2. Add 5.00 g inulin into 0.1 M HAc-NaAc (pH 4.6) buffer, and add the distilled water to make the final volume as 100 mL.

2. Making the standard curve of fructose by DNS method.

Take 6 tubes (with 25 mL volume), and preparing the glucose solution as the below list. And then make these tube went to the DNS reaction, and then the OD₅₂₀ of every tube was recorded. The standard curve of glucose was depicted by using glucose amount as the X-axis, and the OD₅₂₀ value as the Y-axis.

According to the OD₅₂₀ value, the regression equation $Y=aX+b$ was finally got.

| Tube number | 1 | 2 | 3 | 4 | 5 | 6 |
|--|---|-----|-----|-----|-----|-----|
| Fructose solution (mL) | 0 | 0.4 | 0.8 | 1.2 | 1.6 | 2.0 |
| HAc-NaAc buffer (mL) | 2 | 1.6 | 1.2 | 0.8 | 0.4 | 0 |
| DNS solution (mL) | 2 | 2 | 2 | 2 | 2 | 2 |
| Put into the boiling water for exact 5 min, cooled with water, then add water into the tube to make the final volume 25 mL | | | | | | |
| OD ₅₂₀ | | | | | | |

3. Determination of the inulinase activity.

3.1. Sample preparations.

For the powder inulinase formulation, we recommended you dissolved 1 g powder into 100 mL into the 0.1 M HAc-NaAc buffer buffer, stir to make the powder solved.

For the liquid formulation, you can take 1 mL liquid into 0.1 M HAc-NaAc buffer.

3.2. The steps for inulinase activity checking was following:

Take 100 uL diluted inulinase solution and added into 900 uL 5% inulin solution, mixed and incubated in 50 °C water bath for exact 5 min.

Take out the reaction mixture and immediately put into the boiling water for 5 min to inactive the enzyme.

Aliquot 50 uL reaction mixture, and added into 2.0 mL DNS+1.95 mL 0.1 M HAc-NaAc, mixed, and then put into the boiling water for exact 5 min to make the DNS reacted with the reductive sugar. We can see that the color of the liquid becoming dark brown.

Cool the tube with water, and then add the distilled water into the tube to make the final volume was 25 mL.

Then the OD value was determined at OD₅₂₀ nM.

Then referred the standard fructose curve and the regressive equation, the amount of the reducing sugar was calculated ($A = \text{mg/L}$).

3.3. The control experiment.

The inulinase supernatant was firstly put into the boiling water for 5 min to inactive the inulinase.

Described as above, 100 uL inactive enzyme solution was put into 900 uL 5% inulin solution, mixed and then put into 50 °C water bath to react for exact 5 min.

Take out the tube and then put into the boiling water for exact 5 min to inactive the reaction again.

Aliquot out 50 uL reaction mixture, add 2.0 mL DNS+1.95 mL 0.1 M HAc-NaAc, mix together, and then put into the boiling water for 5 min to make the reducing sugar to react with DNS.

Cool the tube by the water to terminate the reaction, and then add distilled water into the tube to make the final volume as 25 mL.

As a control, the solution was taken out as the blank for OD₅₂₀ determination.

Enzyme Unit defined:

One unit was defined as the amount of enzyme need to released one micro mole reduced sugar from the inulin in one minute. The activity of our products was described as the enzyme unite per gram powder or milliliter liquid.

The equation was listed as:

$$U = \frac{A \times 1000}{180 \times 5} \times N$$

U: inulinase activity;

A: The concentration of the fructose in the reaction mixture (mg/mL), This value comes from the regression equation $Y=aX+b$, $A= X$;

180: Molecular weight of the fructose;

5: Reaction time (5 min in this test)

N: Total dilution fold.