

产品名称: **AZD3463**

产品别名: **ALK/IGF1R inhibitor; AZD-3463**

生物活性:				
Description	AZD-3463 is an ALK/IGF1R inhibitor which overcomes multiple mechanisms of acquired resistance to crizotinib. IC50 Value: Target: ALK/IGF1R			
Solvent&Solubility	<i>In Vitro:</i> DMSO : 20 mg/mL (44.55 mM; Need ultrasonic) H₂O : < 0.1 mg/mL (insoluble)			
	Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg
			10 mg	
		1 mM	2.2274 mL	11.1371 mL
		5 mM	0.4455 mL	2.2274 mL
		10 mM	0.2227 mL	1.1137 mL
	*请根据产品在不同溶剂中的溶解度选择合适的溶剂配制储备液; 一旦配成溶液, 请分装保存, 避免反复冻融造成的产品失效。 储备液的保存方式和期限 -80°C, 6 months; -20°C, 1 month。 -80°C 储存时, 请在 6 个月内使用, -20°C 储存时, 请在 1 个月内使用。 <i>In Vivo:</i> 请根据您的实验动物和给药方式选择适当的溶解方案。以下溶解方案都请先按照 In Vitro 方式配制澄清的储备液, 再依次添加助溶剂: ——为保证实验结果的可靠性, 澄清的储备液可以根据储存条件, 适当保存; 体内实验的工作液, 建议您现用现配, 当天使用; 以下溶剂前显示的百分比是指该溶剂在您配制终溶液中的体积占比; 如在配制过程中出现沉淀、析出现象, 可以通过加热和/或超声的方式助溶 1.请依序添加每种溶剂: 10% DMSO→40% PEG300 →5% Tween-80 → 45% saline Solubility: ≥ 2 mg/mL (4.45 mM); Clear solution 此方案可获得 ≥ 2 mg/mL (4.45 mM, 饱和度未知) 的澄清溶液。 以 1 mL 工作液为例, 取 100 μL 20.0 mg/mL 的澄清 DMSO 储备液加到 400 μL PEG300 中, 混合均匀; 向上述体系中加入 50 μL Tween-80, 混合均匀; 然后继续加入 450 μL 生理盐水定容至 1 mL。			