



上海源叶生物科技有限公司
Shanghai yuanye Bio-Technology Co., Ltd
电话: 021-61312973 传真: 021-55068248
网址: www.shyuanye.com
邮箱: shyysw@sina.com

产品名称: **Bay 41-4109 (less active enantiomer)**
产品别名: **Bayer 41-4109 less active enantiomer**

生物活性:

Description	Bay 41-4109 less active enantiomer shows less activity than Bay 41-4109. BAY 41-4109 is a potent inhibitor of human hepatitis B virus (HBV) with an IC ₅₀ of 53 nM.				
IC ₅₀ & Target	IC50: 53 nM (HBV, Bay 41-4109)[1]				
In Vitro	BAY 41-4109 is able to both accelerate and misdirect capsid assembly in vitro. Preformed capsids are stabilized by BAY 41-4109, up to a ratio of one inhibitor molecule per two dimers[2]. BAY 41-4109 is equally effective at inhibiting HBV DNA release and the cytoplasmic HBcAg level, with IC50s of 32.6 and 132 nM in HepG2.2.15 cells, respectively. HBV DNA and HBcAg are inhibited in a dose-dependent manner, indicating that the anti-HBV mechanisms are associated with and dependent on the rate of HBcAg inhibition[3].				
In Vivo	BAY 41-4109 reduces viral DNA in the liver and in the plasma dose-dependently with efficacy comparable to 3TC. BAY 41-4109 reduces hepatitis B virus core antigen (HBcAg) in livers of HBV-transgenic mice. Pharmacokinetic studies in mice have shown rapid absorption, a bioavailability of 30% and dose-proportional plasma concentrations, about 60% in rats and dogs[1]. BAY41-4109 inhibits virus production in vivo by a mechanism that targets the viral capsid[2].				
Solvent&Solubility	In Vitro: DMSO : ≥ 37 mg/mL (93.49 mM) * "≥" means soluble, but saturation unknown.				
	Preparing Stock Solutions	Solvent / Mass / Concentration	1 mg	5 mg	10 mg
		1 mM	2.5268 mL	12.6339 mL	25.2678 mL
		5 mM	0.5054 mL	2.5268 mL	5.0536 mL
		10 mM	0.2527 mL	1.2634 mL	2.5268 mL
	*请根据产品在不同溶剂中的溶解度选择合适的溶剂配制储备液。一旦配成溶液，请分装保存，避免反复冻融造成的产品失效。 储备液的保存方式和期限 -80℃, 6 months; -20℃, 1 month。-80℃ 储存时，请在 6 个月内使用，-20℃ 储存时，请在 1 个月内使用。				
References	[1]. Weber O, et al. Inhibition of human hepatitis B virus (HBV) by a novel non-nucleosidic compound in a transgenic mouse model. Antiviral Res. 2002 May;54(2):69-78. [2]. Stray SJ, et al. BAY 41-4109 has multiple effects on Hepatitis B virus capsid assembly. J Mol Recognit. 2006 Nov-Dec;19(6):542-8. [3]. Wu GY, et al. Inhibition of hepatitis B virus replication by Bay 41-4109 and its association with nucleocapsid disassembly. J Chemother. 2008 Aug;20(4):458-67.				
实验参考:					
	Cellular metabolism is evaluated by MTT colorimetry. HepG2.2.15 cells are plated at a density of 2×10 ³ cells per well in 96-well plates. After 8 d of treatment with different concentrations of each				



上海源叶生物科技有限公司
Shanghai yuanye Bio-Technology Co., Ltd
电话: 021-61312973 传真: 021-55068248
网址: www.shyuanye.com
邮箱: shyysw@sina.com

Cell Assay	antiviral compound, 20 μ L of MTT solution (5 g/L) are added to each well and incubated at 37°C for 4 h. Next, 150 μ L of DMSO is added and stirred for 10 min to dissolve the crystals. Absorbance values are recorded at 490 nm by using an ELISA reader. The MTT values are calculated using the curve regression equation ^[3] .
Animal Administration	Mice: The HBV-transgenic mice are used in the study. Compounds (BAY 41-4109) are formulated as a suspension in 0.5% Tylose and administered per os to mice two times/day for a 28 day period. The 0.5% Tylose serves as a placebo. Six hours after the last treatment, the animals are sacrificed and livers are removed and immediately frozen for subsequent analysis. Blood is obtained by cardiac puncture of the anesthetized animals ^[1] .
References	[1]. Weber O, et al. Inhibition of human hepatitis B virus (HBV) by a novel non-nucleosidic compound in a transgenic mouse model. <i>Antiviral Res.</i> 2002 May;54(2):69-78. [2]. Stray SJ, et al. BAY 41-4109 has multiple effects on Hepatitis B virus capsid assembly. <i>J Mol Recognit.</i> 2006 Nov-Dec;19(6):542-8. [3]. Wu GY, et al. Inhibition of hepatitis B virus replication by Bay 41-4109 and its association with nucleocapsid disassembly. <i>J Chemother.</i> 2008 Aug;20(4):458-67.

源叶生物